

## 8 Key Aspects of Avazzia Micro-Current Signal

1. Signal triggers brain's repair response.
2. Increased perfusion at site.
3. Neuropeptides released and production ramped up. Results in pain relief and speeds healing.
4. Reduces Inflammation, swelling, edema - rapidly.
5. Increases ATP in Mitochondria.
6. Stimulates C-fibers- insures continued healing and pain reduction.
7. Increases Lymphatic drainage.
8. Continued use is Accumulative- pain relief lasts longer between uses; rate of healing maximized.

How & Why does the use of Micro-current help someone look younger?

1. As a result of increased perfusion, the skin experiences a renewed and rejuvenated look.
2. As a result of the underlying muscle stimulation, the skin will take on the appearance of greatly increased smoothness- EASILY noticeable. .
3. The signal activates acupuncture points on the face that affect the entire body and generates a feeling of relaxation.
4. Micro-current causes an increase in energy, which may be noticed as having a little more "spring in the step".
5. Repeated use provides **accumulative** benefits. The results look better and last longer. One single use lasts up to 7 to 10 days. Twice weekly use maintains peak performance and appearance.
6. Results are immediate
7. Production of collagen and realignment of collagen in skin helps to reduce the appearance of scars/scar tissue.
8. Technology is non-toxic, non-invasive, has no side-effects, simple to use, can even be used over make-up, is portable (hand-held), can be used anywhere and anyplace user sees fit.

